## Celtra® Duo

Correction Porcelain, Universal Stain and Overglaze

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Mode d'emploi

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**Indications:**

- Celtra Duo Correction porcelain is indicated for minor repairs of Celtra Duo restorations only, including adjustments of contacts and occlusal surfaces.
- Dentsply Sirona Universal Stain and Overglaze are indicated for staining and glazing Celtra Duo restorations. They can also be used for staining and glazing with most lithium silicate based restorations, zirconia based restorations and most of the Dentsply Sirona Porcelains.

**Contraindications:**

Only the indications listed above are suitable. Porcelains other than Celtra Duo Correction Porcelain may not have the correct coefficient of thermal expansion, and therefore should not be used on Celtra Duo Blocks.

**Warnings/Adverse Reactions:**

Consult Dentsply Sirona Universal Stains and Overglaze Safety Data Sheets (SDS) for safe handling and usage. Use only according to the indications listed above.

**Precautions:**

Refer to the special notes throughout the instructions below for best practices for using Dentsply Sirona Universal Stains and Overglaze.


- This product is intended to be used only as specifically outlined in these Directions for Use. Any use of this product inconsistent with the Directions for Use is at the discretion and is the sole responsibility of the practitioner.

**Storage conditions:**

Store materials in packages with the lid tightly sealed. Do not leave open bottles where dust and debris may contaminate the porcelain. Replace cap when not in use. Avoid vibration and heat. Keep dry. Shake well before use. The lot number is printed on each bottle; please refer to the lot number in all correspondence. Do not use after expiration date.

**Crown Preparation:**

- Celtra Duo restorations require adequate preparation reduction and restoration thickness. Insufficient wall

thickness may lead to premature failure.

\*1.5-2 mm for frameworks that are exclusively polished

- Celtra Duo restorations are designed to be polished and/or glazed before insertion. Direct insertion without polishing or glazing may lead to excessive wear.
- Use only in well ventilated areas.

**Finishing the restoration**

Without the optional staining and/or glaze firing

1. Clean the restoration in an ultrasonic water bath or with a steam cleaner. Carefully blow-dry the restoration in an air stream.
2. Try in the final restoration to examine its fit. If necessary, adjust the proximal and occlusal contacts.
3. When making adjustments to the ceramic material, make sure to avoid overheating (e.g. by providing irrigation). Ideally, all adjustments should be made with water-cooled instruments.
4. Diamond polishing bodies (< 60µm) are recommended for polishing the occlusal surfaces. The use of tungsten carbide burs is not recommended.

5. Smooth the remaining surfaces of the restoration with a rubber polisher at low speed and minimum pressure.
6. Thoroughly clean the restoration.
7. If milled restorations require no corrections, proceed to optional staining and glazing instructions.
8. If milled restorations require corrections, proceed to Application of Celtra Duo Correction Porcelain instructions.

**Instructions for Application of Celtra Duo Correction Porcelain to Celtra Duo CAD/CAM restorations**

**Procedure:** Mix Celtra Duo Correction Porcelain with distilled water or Dentsply Sirona Modeling Liquid U to desired consistency. Apply porcelain where applicable. Fire to complete the restoration following the procedures and firing programs listed below:

**First Glaze & Celtra Duo Correction Porcelain on a Crown** - Apply Dentsply Sirona Universal Overglaze, over the entire Celtra crown then apply Celtra Duo Correction porcelain to the required areas and fire at 820°C according to 1st correction firing program a) in Table 1. Polish to a high gloss, if needed.

**Second Glaze:** If after first firing the required glaze appearance is not achieved, re-apply glaze, and fire according to 2nd & subsequent glaze firing program in Table 1.

**Table 1. Correction/Glaze Firing**

	a) 1st correction firing/ b) 2nd correction firing	2nd & subsequent glaze firing
Pre-drying	2:00 min or 4:00 min when SuperPeg II™ used	2:00 min or 4:00 min when SuperPeg II™ used
Drying <sup>1</sup>	2:00 min	2:00 min
Preheating	2:00 min	2:00 min
Start Temp	500°C	500°C
Heating rate	55°C/min	55°C/min
Final Temp	a) 820°C b) 810 °C	770°C
Vacuum	on 1:00 min	off
Hold time	0:30 min	1:30 min
Cooling <sup>2</sup>	3:00 min	3:00 min

<sup>1</sup> Drying time depends on the type of furnace

<sup>2</sup> This cooling period of three minutes is only required for objects placed on firing pins or when using an auxiliary firing paste. If the object is placed directly on firing cotton with a firing support, no cooling phase is required.

**Second Application of Glaze & Celtra Duo Correction Porcelain (or application of Celtra Duo Correction Porcelain after previous glazing):** If a second application is needed: apply Dentsply Sirona Universal Overglaze, over the entire Celtra Duo crown, then apply small amount of Celtra Duo Correction Porcelain to the desired areas. Fire at 810°C according to 2nd correction firing program b) in Table 1.

Note: Never fire the 2nd correction with 820 °C. Contour the areas where Celtra Duo Correction Porcelain was applied, if needed, and polish to a glaze finish.

**Third Glaze (optional):** If after firing the required glaze appearance is not achieved, re-apply glaze over the entire unit, and fire following 2nd & subsequent glaze firing program in Table 1.

### Instructions for optional staining and glaze-firing of Celtra Duo CAD/CAM restorations

Custom staining and glazing of Celtra Duo restorations can be performed using Dentsply Sirona Universal Stains and Overglaze instructions below.

- For best results, any sanding marks on the surface should be removed with a suitable ceramic-coated rubber polishing wheels or polished smooth at moderate speed (not exceeding 8,000 rpm). The setting of the furnace parameters depends on the type of furnace, so the programming may vary slightly.

### Procedure: For Dentsply Sirona Universal Stains and Overglaze application the following steps must be observed:

- Before applying the stains and glaze, the restoration must be clean and free of grease. Clean the surface of the restoration with a steam cleaner or submerge for 10 minutes in an ultrasonic cleaner with distilled water. Any contamination after cleaning must be prevented.
- Place the cleaned crown on the prepared DENTSPLY Prosthetics die material. To achieve the correct shade, prepare the die material according to the shade recommendations in a Table 2 below.

**Table 2. Recommended Die Material Selection**

Shade	A1	A2	A3	A3.5	B2
Die Material	F1	F12	F10	F9	F11

- Mix the mass thoroughly with spatula. Do not use a metal spatula. Withdraw desired amount of Dentsply Sirona glaze and place it on the mixing palette. For a thinner consistency dilute the material with Dentsply Sirona Universal Stain and Overglaze Liquid.
- Apply a thin layer of glaze to the entire crown surface using a brush. Apply a thin layer of stain in the gingival area and spread with a brush toward the incisal region (or as needed). Check the shade using the shade tab. Use the brush to increase or decrease the stain amount to achieve the desired shade.
- Again using the brush, apply a thin layer of incisal stain to the incisal edge. Check the shade using the shade ring. Use the brush to increase or decrease the stain amount to achieve the desired esthetic incisal effects.

#### Notes:

- A more intensive shade effect can be achieved by repeating cycles of applying and firing the material. However, the use of too many layers of stain may result in an unnatural appearance.
- The cusps and fissures can be individually characterized with stains.
- The basic shade is determined based on the Vita shade groups (A, B, C, D) (see Table 3).

**Table 3. Celtra Duo Body & Enamel Stains Recommendations for Celtra LT Blocks**

Shade	Block	Body Stain	Incisal Stain
A1	A1	A wash of Stain 1, or just glaze	i1
A2	A2	A wash of Stain 1, or just glaze	i1
A3	A3	A wash of Stain 1, or just glaze	i1
A3.5	A3.5	A wash of Stain 2, or just glaze	i2
A4	A3.5	<b>1st firing:</b> Stain i1 wash over the entire crown, including incisal area <b>2nd firing:</b> <b>Body:</b> 50/50 mix of Sunset & Stain 3 <b>Incisal:</b> Stain i2	i2
B1	A1	Stain 2	i1
B2	B2	A wash of Stain 2	i1
B3	B2	1 firing of Stain 2	i1
B4	B2	2 firings of Stain 2	i1
C1	B2	A wash of Stain 3	i2
C2	B2	Stain 3	i2
C3	B2	<b>1st firing:</b> 2/3rd Stain 3 & 1/3rd Stain 2 mixed together and wash over the entire crown, including incisal area <b>2nd firing:</b> <b>Body:</b> 50/50 mix from 1st fire & glaze <b>Incisal:</b> 50/50 mix of i1 & glaze	50/50 mix of i1 & glaze
C4	A3.5	<b>1st firing:</b> 50% Stain 3 & 50% i1 mixed together and wash over the entire crown, including incisal area <b>2nd firing:</b> <b>Body:</b> 50/50 mix from 1st fire & glaze <b>Incisal:</b> 50/50 mix of i1 & glaze	50/50 mix of i1 & glaze
D2	A1	Stain 4	i1
D3	A2	Stain 4	i1
D4	B2	<b>1st firing:</b> 2/3rd Stain 3 & 1/3rd Stain 2 <b>2nd firing:</b> Olive	i2


### Use of glazes and stains by third-party manufacturers

Only the native range of glaze and stains (Dentsply Sirona Universal Stain and Glaze) in conjunction with the appropriate liquid is recommended for customizing and glazing Celtra Duo. Glazes and stains by third-party manufacturers may require other firing temperature and cycles than those recommended for Celtra Duo.

### Firing—Dentsply Sirona Universal Stains and Overglaze Only

**General:** Firing objects must be placed on a carrier within the furnace. **Placing the objects on plain firing pad with firing tray or on firing pad with a platinum foil is recommended.**

Alternatively, scale-free rounded metal pins (e.g. platinum wire) or firing pins may be used, or firing cotton for direct support or in combination with thin metal or ceramic pins. Make sure that the pins do not touch the restoration. If you prefer using an auxiliary firing paste, Harvest Dental SuperPeg Putty II™ refractory paste is recommended (available separately, see manufacturer's complete Directions for Use), applied loosely to the inside of the crown and never beyond the crown margin. Thick ceramic pins or third party auxiliary firing pastes may result in cracks and are therefore not recommended for use in connection with Celtra Duo restorations.

When using firing pins, make sure that no sharp pin ends touch the restoration. Firing paste should always be applied in small quantities and should not completely fill the inside of the restoration.

- Always use a firing pad for posterior units and/or looped Dentsply Sirona All-Ceramic firing pins with a small amount of Super Peg Putty, placed onto a honeycomb sagger tray. Do not use other types of firing pins.

#### Firing paste

- If a firing paste is preferred, then only SuperPeg II™ is recommended.
- Other firing pastes may lead to restoration cracking.

### General firing recommendations

- Follow glaze/stains firing programs listed in Table 4.
- The optional glaze firing will increase the flexural strength of Celtra Duo to 370 MPa. To increase the flexural strength it is only necessary to fire the restoration up to 820°C. The application of glaze is an option but not necessary. In case of firing without glaze skip the Pre-drying and drying steps of the firing program. Make sure to polish the restoration first before firing.
- Starting temperature 500°C
- Place the object on the firing tray/firing table of the furnace.

**Table 4. Glazing/Stains Firing**

	1st glaze/stain firing	2nd & subsequent glaze/stain firing
Pre-drying	2:00 min or 4:00 min when SuperPeg II™ used	2:00 min or 4:00 min when SuperPeg II™ used
Drying <sup>1</sup>	2:00 min	2:00 min
Preheating	2:00 min	2:00 min
Start Temp	500°C	500°C
Heating rate	55°C/min	55°C/min
Final Temp	820°C	770°C
Vacuum	off	off
Hold time	1:30 min	1:30 min
Cooling <sup>2</sup>	3:00 min	3:00 min

<sup>1</sup> Drying time depends on the type of furnace

<sup>2</sup> This cooling period of three minutes is only required for objects placed on firing pins or when using an auxiliary firing paste. If the object is placed directly on firing cotton with a firing support, no cooling phase is required.

- Additional glaze firings may be performed at 770°C in order to accentuate the shade or to correct it with glaze, or increase the gloss. It is necessary to coat the entire surface with Dentsply Sirona universal glaze.
- Finally, check the shade against a shade tab and adjust if necessary.
- Note: After the first firing, the restoration may not appear sufficiently glazed. In this case, apply a thin layer of Dentsply Sirona glaze to the entire surface with a brush. After applying the glaze, perform a second stain/glaze firing as defined in the firing recommendations. Additional stain/glaze firing cycles can be performed using the same parameters.

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Super Peg II™ is not a registered trademark of DENTSPLY International, Inc.
Super Peg II™ n'est pas une marque déposée de DENTSPLY International, Inc.
Super Peg II™ non è un marchio registrato di DENTSPLY International, Inc.
SuperPeg II™ ist keine eingetragene Marke von DENTSPLY International, Inc.
Super Peg II™ no es una marca comercial registrada de DENTSPLY International, Inc.
DENTSPLY International Inc.